

# COMON



**safety**

MOVING ON THE BATHROOM

**help**

WALKING AROUND

**move support**

ON BED

**enjoy**


DAILY LIFE

**share**

FREELY EDLER HEATHCARE

**catalogue**<sup>2025</sup>







## ABOUT **comon**

---

Comon is...

located in Zhejiang province, Founded in 2012. Employs more than 100 people, sophisticated modern production equipment, Specializing in the production of wheelchair, walker, walking sticks, Handy Grabber, walking aid and other medical equipment.

it is a company one-stop service of development, production and sales. Products are widely used in medical rehabilitation, disability and environmental protection industries. Our company has paid attention to technological progress and scientific management for a long time. Furthermore, we actively learn from domestic and foreign advanced technology, management experience, and product innovation, to obtain a large and steady market share in fierce competition.

For many years, enterprises have been committed to the development of new products to meet the growing needs of the market. We offer a reliable product quality, and win customers' good reputation. Our products exported to Europe, the United States, Japan, South Korea, South-east Asia as well as Australia and other markets. Company have been certified the BSCI, SEDEX and ISO9001.2015, and also have obtained many internationally recognized certificates for all products, which are complied with relevant European and American approval standards such as CE, MDR, FDA.

We believe "there is no best, only better". Enterprise has always implemented the principles of "customer satisfaction is our eternal pursuit", holding pragmatic cooperation attitude and flexible business strategy and dedication with the new and old customers sincere cooperation.



5



31



43

19





# TABLE OF **content**

---

P.01-02 **About Comon**

P.03-04 **Table of content**

P.05-18 **Bathroom & Toilet Safety**

P.19-30 **Mobility & Travel**

P.31-42 **Sleeping & Living Aid**

P.43-50 **Accessories & Printing  
design**



## **More about Comon**

Go to:  
[www.cncomon.com](http://www.cncomon.com)

## **Comon medical appliances Co.,Ltd**

No.41 Develop Road Wuyi Zhejiang, China  
+86.15088231753  
[sales@cncomon.com](mailto:sales@cncomon.com)







# bathroom & toilet safety

---

Bathrooms can be a dangerous place, especially for older adults and children. Slippery floors, sharp edges, and hot water are just a few of the hazards that can lead to accidents and injuries.

In this article, we will provide you with a comprehensive guide to bathroom safety, including tips for preventing accidents, choosing the right equipment, and creating a safer bathroom environment for everyone.

# bath & shower seat

A shower chair allows a person to sit in the shower while bathing. They are plastic, have at least four nonslip feet, and generally fit fully within the tub or shower. Chairs are available with or without a back and are good for those who may struggle with balancing in the slippery environment.



## shower stool

This is an ordinary shower chair, made of lightweight aluminum material and equipped with anti slip foot covers. The backrest style makes you more comfortable

1. SHOWER STOOL  
2001

2. SHOWER STOOL BACKREST  
2004



## ergonomic seat shower stool

The ergonomic seat offer more convenient cleaning of the body

1. ERGONOMIC SHOWER STOOL  
2007

2. ERGONOMICS HOWER STOOL BACKREST  
2008



## small size shower stool

Small size design provides greater possibilities for e-commerce

1. SMALL SIZE SHOWER STOOL BACKREST  
2009

2. SHOWER STOOL BACKREST  
2010





## handrail shower stool

The armrest design makes it easier for elderly people to get up and sit down, and the rubber pad prevents the body from slipping off the stool surface. Strengthening pipe design to improve product stability and load-bearing capacity.

### 1. SHOWER STOOL WITHOUT HANDLE

2012

### 2. SHOWER STOOL WITH HANDLE

2011

### 3. STRENGTHEN SHOWER STOOL WITH HANDLE

2028



## round seat shower stool

---

The circular stool seat can sit in any direction

1. ROUND SEAT SHOWER STOOL  
2015



## detachable shower stool

---

Detachable legs can provide more convenient and fast assembly without affecting packaging

1. DETACHABLE SHOWER STOOL  
2029



## rotating shower stool

---

Rotate 360 degrees without moving

1. TOTATING SHOWER STOOL  
2022



## **round shower seat**

The circular stool surface design has no restrictions on position and direction, providing greater possibilities for narrow spaces and bathtubs. 13.4"/20.9" height adjustment, also made of lightweight material and anti slip foot covers

## step stool

As you age, it becomes more difficult to raise your legs and feet. The risk of falling when these people enter the bathtub comes with it. The standing stool is a great solution and provides a range of adjustable sizes for selection

1. **STEP STOOL**  
2016

2. **STEP STOOL ADJUSTABLE**  
2017





## short grip bathtub seat

---

The bathtub seat is a type of stool that spans over a bathtub. Small in size and easy to install, it is loved by many bathtub enthusiasts. You can sit on it to take a shower, or have a glass of red wine or food while soaking in the bathtub

1. SHORT GRIP BATHTUB SEAT  
2013



## angeld bathtub seat

---

This angled design not only has the normal function of a bathtub seat, but can also be inverted and assembled into a standing stool. It can be adjusted in two sizes, 28.3 inches and 32.3 inches, depending on the width of the bathtub. Extending the handle cover allows for a larger grip on the ground

1. ANGELD BATHTUB SEAT  
2014







## detachable transfer bathtub seat

---



The Plastic Tub Transfer Bench with Backrest has four individually-adjustable legs to sit level at the optimal height for each individual regardless of the height difference of the tub floor and bathroom floor. The bench provides a stable, secure surface to assist individuals with mobility limitations in and out of the tub, allowing them to retain more of their freedom and independence. The back can be attached to either side of the bench, and both arms are independently removable, allowing the chair to be left- or right-handed to suit any bathroom orientation.

The most important thing is that this product adopts a fast assembly design, which solves the problem of difficult assembly of this product. With just a few pipes assembled, this product can be quickly completed

1. TRANSFER BATHTUB SEAT  
2019

# transfer bathtub bench

---

The three piece shower bench is an extended shower bench that is mainly convenient for people with insufficient lower limb strength to support themselves in and out of the bathtub. One end of the product is placed in the bathtub, and the other end is on the ground outside the bathtub. It has a large suction cup design that firmly sucks the ground, allowing users to quickly enter and exit the bathtub through the armrest



1

## 1. TRANSFER BATHTUB BENCH

2018

## 2. BATH BENCH SHOWER BRUSH

201801



2



# toilet chair & handrail

The biggest problem for elderly, pregnant women, and some people with limited mobility is the inability to sit down and get up on the toilet. They must rely on a raised toilet to sit down safely and stand up on the handrail next to the toilet after finishing. This device effectively solves the troubles of these people



# toilet seat & chair

---

The toilet chair adopts a quick assembly design, with left and right 2-speed adjustment suitable for various sizes of toilets, and 8-speed height adjustment suitable for people of various heights. The 300 pound load-bearing design meets the needs of most people. The seat cushion is made of very soft material and offers backrest styles for selection

- 1. TOILET CHAIR  
2030
- 2. TOILET SEAT  
2031



# toilet handrail



1

1. TOILET HANDRAIL  
2032

2. TOILET HANDRAIL BAG  
203201



2



# toilet handrail

---

Simple toilet armrest

1. TOILET HANDRAIL  
2026



# steel toilet handrail

---

The iron material makes this product more stable and load-bearing, and the clamping devices on both sides can firmly clamp the armrest onto the toilet without shaking. The width and height can be adjusted as needed. The part that fits the toilet adopts a curved design to make the product fit better

1. STEEL TOILET HANDRAIL  
2033



# toilet handrail

---

Simple toilet armrest

1. TOILET HANDRAIL  
2025



# steel toilet handrail

---

The steel toilet armrest is welded with round pipes, and the two curved pipes are more firmly attached to the toilet. Blue and white designs are fashionable and beautiful

1. STEEL TOILET HANDRAIL  
2034







# mobility & travel

---

As individuals age, many develop chronic, complex illnesses, including deconditioning, and may need to use a mobility assistive device. Assistive devices such as canes, crutches, walkers, and wheelchairs can help to alleviate the effects of mobility limitations, providing improved independence.

According to data from the National Health and Aging Trends Study, 29.4% of adults 65 years and older reported using assistive devices within the previous month when outside the home, and 26.2% reported using them inside the home. An adult should understand how to select the appropriate assistive device

## quad cane

Aquad cane has four small feet extending from the base that make it more stable. It can handle more weight bearing and increases the base of support and stability. This design, which can be purchased with a small or a large base, is an advantage when navigating stairs. Also, when you sit or need to free your hands, the cane stands by itself so you don't have to bend down to pick it up. This kind of cane is especially helpful for seniors and anyone affected by vertigo, partial paralysis, hemiplegia, stroke, or who are recovering from hip or knee replacement surgery.



### 1. QUAD CANE SMALL BASE

1006

### 2. QUAD CANE LARGE BASE

1005

## offset cane

An offset cane handle has a more ergonomic grip than a standard crook handle. Via the bend in the cane shaft, these handles are positioned over the length of the shaft, allowing for even weight distribution along the shaft. With the weight centered, pressure is taken off the hand and wrist. These padded handles reduce strain on your wrist and provide a comfortable grip, coming in multiple styles and shapes





## decarative offset cane

The curvature design of this cane determines its difficulty in printing, We have an exclusive transfer printing method to print on this offset cane. Can achieve no color fading within 2 years and provide customized decors. This cane can also be made in solid colors, such as the popular glossy black, blue, and pink colors to choose from.

1. **OFFSET CANE 515**  
1003515

3. **OFFSET CANE 516**  
1003516

5. **OFFSET CANE BLU**  
1003BLU

7. **OFFSET CANE PNK**  
1003PNK


2. **OFFSET CANE 517**  
1003517

4. **OFFSET CANE 518**  
1003518

6. **OFFSET CANE 514**  
10035147

8. **OFFSET CANE 564**  
1003564





## circular hook cane & seat cane



A tourist handle is also known as a crook or a “J” handle cane, a shepherd’s cane, or a hospital cane. These canes are commonly used to help alleviate general pain from standing and walking. They also are counted on to assist with balance. They are simple and traditional and very common.

## walking seat cane

For anyone who enjoys venturing out and wants to accommodate the need to stop and rest, a walking cane with a seat is a game changer. Also called a walking cane chair or a seat cane, this type of cane greatly contributes to mobility and independence. The seat eliminates barriers by making longer outings possible and comfortable for cane users.



1. WALKING SEAT CANE  
10121018



## crook walking cane

1. CROOK WALKING CANE WITH PVC HANDLE  
1004SL
2. CROOK WALKING CANE WITH EVA HANDLE  
1004BRN

This is a very traditional product, with the usual popular colors being silver and brown .The handle is available in PVC and EVA for selection



## soft-grip folding cane

---

1. **SOFT GRIP FOLDING CANE PURPLE**  
1001PUR
2. **SOFT GRIP FOLDING CANE PRINTING**  
1001DBU
3. **SOFT GRIP FOLDING CANE PRINTING**  
1001LPD
4. **SOFT GRIP FOLDING CANE BLACK**  
1001BLK

The soft rubber handle design makes this product more comfortable to use, and the surface sandblasting process makes the product more textured. Compared to most men's solid color choices, women are more willing to choose our printing process. Printing technology makes products more fashionable.

For the convenience of storage and carrying, we have designed a self-pulling cloth bag for choice.



A person wearing a white t-shirt and white pants is walking on a paved path. They are holding a black walking cane with a Fritz handle. The cane has a black strap attached to the handle. The person is also wearing white sneakers. In the background, there are green plants and a stone wall.

## **fritz handle walking cane**

A Fritz handle can be used on a variety of canes, including single tip canes, folding canes, and those made of just about any material. They resemble the letter T when viewed from the side, and were designed by a German Count in the 1700s to make using a cane more comfortable for someone with arthritis. This open-ended styled handle allows for more room for the user's fingers, and the slight curve rests comfortably in the palm. Sometimes this kind of cane is called a crutch handle cane. There are currently two main materials available: wood and plastic



## full soft grip folding cane

The derby grip handle features innovative technology that decreases shock and reduces fatigue in the hand and upper body. The ergonomic cane handle is adaptable for right- & left-handed users.

1. **FULL SOFT GRIP FOLDING CANE PRINTING**  
1015LOTUS
2. **FULL SOFT GRIP FOLDING CANE BLACK**  
1015BLK
3. **FULL SOFT GRIP FOLDING CANE PRINTING**  
1015RED VEIN
4. **FULL SOFT GRIP FOLDING CANE PRINTING**  
1015GRN VEIN
5. **FULL SOFT GRIP FOLDING CANE PRINTING**  
1015PUR VEIN





## wooden handle folding cane


The wooden handle folding cane, with a glossy black oxidation, is currently the most widely used and widely used cane in the market. It is suitable for all people in need, can be used to assist walking, and can also be used to decorate oneself. Folded up for easy portability, the aluminum alloy material makes it only weigh about 350g.

### 1. WOODEN HANDLE FOLDING CANE 1016



**STORAGE BAG**  
101603

When no crutches are needed, the cloth bag can store the crutches for easy carrying and protect them from damage when dropped



**WRIST LOOP**  
101601/101602

The wrist strap can be fastened on the hand to prevent the cane from slipping or falling off



**ICE TIP**  
101604

Crutch accessories, convenient for walking on ice or slippery muddy roads

# crutch

Injuries can happen to anyone at any time. If you go to a hospital for treatment, the doctor will probably plaster the fractured area of your bone. After that, you will probably need a pair of crutches. When you injure the lower limbs of your body, such as your ankle or knee or any other joint, you might feel uncomfortable while walking. You need forearm and elbow crutches that can be adjusted to reduce the discomfort caused by the plaster. This type of plaster will also help you in walking when you will be recovering from the injury.



## soft grip crutch

---

For someone who is injured in the leg and unable to move normally, all of their strength is focused on the palm of their hand while moving. Therefore, an excellent elbow is particularly important. The palm support part of this elbow is made of soft rubber, which effectively shares the pain caused by weight. The circular design of the arm effectively wraps around the handle, making movement more flexible

### 1. SOFT GRIP CRUTCH 1011



## crutch

---

Aluminum elbow (forearm) with gradual height adjustment. crutch. It is designed to deliver safe support for persons who need additional stabilization during walking. It is available in gray and blue.

### 1. CRUTCH BLUE 1010BLU







# sleeping & living aid

---

According to the statistics 61% of people would like to stay in their own home if they reached a point where they could no longer live safely by themselves.

With the help of independent living products, as well as medical technology, and some assistance from other people, seniors can retain a certain level of independence. By increasing independence, and thus dignity, elderly are likely to feel more confident and happy. Which is why these assistive devices make great gifts too



## 2-piece design bedside table

the hospital bed table also provides a large, stable surface for mealtimes and other activities when recovering from an illness, injury or surgery. Height adjustable, the over bed table utilizes four swivel castors to easily position it for left- or right-handed use at any bed, sofa or chair. Constructed with a durable steel frame with a laminate top, the bedside table is easy to assemble with illustrated step by step instructions.

1. 2-PIECE DESIGN BEDSIDE TABLE  
3001





## bedside table

---

Classic style bedside table, quick tool free assembly

1. **BEDSIDE TABLE WHITE**  
3002WHT
2. **BEDSIDE TABLE GINGER**  
3002GIN



## square base bed table

---

Flexible movement of small base bedside table, suitable for smaller gaps

1. **SQUARE BASE BED TABLE WHITE**  
3003WHT
2. **SQUARE BASE BED TABLE GINGER**  
3003GIN



**CIRCULAR BED HANDLE**

3004

One of the popular styles, with a simple structure. Lightweight use



**BED HANDLE**

3005

Popular e-commerce styles with small and fully functional packaging, and a handle added on the side for the convenience of physically weak people to hold onto and wake up



**STRENGTH BED HANDRAIL**

3006

This product is designed with a larger part that is fixed to the bed, and the handle part has a greater load-bearing capacity




## bed support

Getting up and turning over are simple actions that are not easy for the elderly or patients lying in bed. A product that provides auxiliary support is particularly important. A bedside table is a product that prevents elderly people or patients from falling during night rest and provides support for those who have difficulty moving while lying in bed







**STEEL BED HANDLE**  
3007

---

Traditional European style steel design



**ALUMINIUM BED HANDLE**  
3008

---

Lightweight aluminum handrails and cloth bag design make it convenient to store some commonly used items



**BED HANDLE**  
3009

---

E-commerce style, small packaging volume, suitable product size




## bed backrest

---

When you need to lean on the bed to read or work, or when people with limited mobility need to do some tasks done in bed, maintaining unsupported movements for a long time can cause back pain and lower back pain. the bed backrest can provide consistent comfort and relief

1. **BED BACKREST**  
3010





# urine bottle holder

Getting up to use the toilet at night time can be a struggle for most, particularly bedridden people, sufferers of mobility issues, and elderly individuals. The Wire Urinal Bottle Holder for Bed Headboards, can help to prevent this need hanging neatly over the majority of headboards, for use day and night without risk of spills. Featuring a plastic coating for easy cleaning, this design upholds hygiene at home and in care while promoting convenience for patients and carers

## 1. URINE BOTTLE HOLDER

3011

## 2. URINE BOTTLE HOLDER WITH LID

3012



## sock aid

Bending down to tie shoelaces or pull on a pair of socks gets harder with age; dressing yourself with the Sock Aid helps prevent over-extension from bending down to pull on your favorite pair.





## cloth sock aid

---

Blue Oxford cloth allows feet to smoothly slide into socks.  
White cotton fabric can firmly support socks

### 1. CLOTH SOCK AID 3013



## plastic sock aid

---

PE plastic smooth surface, nylon  
rope and sponge handle

### 1. PLASTIC SOCK AID 3014

## steel sock aid

---

Stable steel structure, able to stand  
on the ground and complete the  
process of putting on socks with  
feet extended in

### 1. STEEL SOCK AID 3015





# handy grabber

---

No shelf is out of reach when you use the Grip 'n Grab. With this arm extension, you can get items down from high spots or from the floor without losing your grip or suffering a tragic fall.



## foldable handy grabber

---

This is a very popular soft rubber design that can firmly grip items of various shapes and is easy to carry when folded. The handle is also a comfortable soft rubber.

1. **FOLDABLE HANDY GRABBER**  
3016



## magnetic reacher

---

The scissor shaped design allows for picking up very small things. The strong magnetic field on the head can easily attract keys and coins.

1. **MAGNETIC REACHER**  
3017



## square handy grabber

Square groove design, compact and easy to carry. The opening can be adjusted according to the object being grabbed. And there are shoe pull styles to choose from, making it a great helper for people with limited mobility



1. **SQUARE HANDY GRABBER FOLDING**  
3019
2. **SQUARE HANDY GRABBER UNFOLDING**  
3018
3. **SQUARE HANDY GRABBER SHOES HORN**  
3020



**HIP KIT**  
3021

Basic daily necessities set,include handy grabber,dressing aid,shoes horn ,scrubbing bar and sock aid



**HIP KIT**  
3022

Basic daily necessities set,include handy grabber,dressing aid,steel shoes horn ,scrubbing bar,sock aid and shoelaces



**HIP KIT**  
3023

Basic daily necessities set,include folding handy grabber,dressing aid,shoes horn ,scrubbing bar and cloth sock aid







# accessories & printing design

---

All products are consumables and wear parts, so we provide replacement parts for all consumables to extend the service life of products and reduce their usage costs, better serving product users

# antiskid rubber tip

Antiskid tips are essential for stability and traction. Rubber tips, also known as a cane ferrule, are almost always on a walking cane, stick and shower stool. Choosing the right size, rubber material, and tread is essential to purchasing a antiskid tip. Ensuring they are not worn out is also fundamental to safely using a cane or chair and avoiding accidents.



## trusty cane tip



Trusty Cane is the one cane you can trust to provide secure, reliable support on any surface. Trusty Cane's triple tread base is designed to improve stability for a more natural stride wherever you go. Ordinary cane tips can slip when you walk, but Trusty Cane's unique pivot lets the cane move with you while the triple tread base holds its grip.

### 1. TRUSTY CANE TIP 4001

## trusty cane tip



Trusty Cane is the one cane you can trust to provide secure, reliable support on any surface. Trusty Cane's triple tread base is designed to improve stability for a more natural stride wherever you go. Ordinary cane tips can slip when you walk, but Trusty Cane's unique pivot lets the cane move with you while the triple tread base holds its grip.

### 2. TRUSTY CANE TIP 4002



## suction cup tip

---

Large suction cup footrest, firmly suction on the ground, used for bathroom shower benches

### 1. SUCTION CUP TIP 4003



## large base tip

---

Large base footrest, commonly used in products such as shower benches, toilet chairs, etc

### 1. LARGE BASE TIP 4004



## trusty cane rubber tip

---

Rubber four legged leg covers are suitable for most crutches and can be replaced at any time. The weight of 400g provides excellent standing support for most crutches

### 1. TRUSTY CANE RUBBER TIP 4005



### WALKING CANE TIP 4006

---

Regular crutch leg covers are available in black, gray, and coffee colors. The rubber wear-resistant material is suitable for various weather and ground conditions



### SHOWER CHAIR TIP 4007

---

The most widely used regular shower bench leg cover is sturdy and durable. Rubber material exhibits excellent grip on wet surfaces



### SHOWER CHAIR TIP 4008

---

Angled shower stool foot cover, suitable for angled shower stools, providing better anti tipping ability for the stool





**SOFT GRIP HANDLE**  
5001

---

Soft rubber handle, currently one of the most commonly used plastic handles, provides excellent palm cushioning with its soft rubber design



**SOFT GRIP HANDLE**  
5002

---

Fully covered adhesive handle, one of the popular styles of cane handle



**PLASTIC HANDLE**  
5003

---

The most conventional plastic handle with a frosted surface design. Good strength and high load-bearing capacity





## left and right handle

The left and right handles are unidirectional crutches designed based on human hands. He is only suitable for left or right hand use, but provides a good fit for the handle, making it tighter and more comfortable to hold in the hand



1. **RIGHR HANDLE**  
500401R  
2. **LEFT HANDLE**  
500401L

3. **LEFT HANDLE**  
500402L  
4. **RIGHT HANDLE**  
500402R

## wood handle

Conventional wooden handle



1. **WOOD HANDLE**  
5006

## wooden handle

The most common and popular wooden handle cane, suitable for all crutches. Wood material makes crutches more upscale



## bird handle


Bird shaped handle


1. **WOODEN HANDLE**  
5005

1. **BIRD HANDLE**  
5007

## walking cane handle

We offer a vast selection of different types of walking canes for our customers, and we want to ensure you are comfortable with your selection. It will help you make the right selection. Some cane handles are designed for comfort while others are made to be stylish and look great. It all depends on what you want from your cane.







**COMON medical appliances LTD**

NO.41 Development Road

Baihuashan Industry Zone

ZHEJIANG CHINA

+86(579)87928238

[sales@cncomon.com](mailto:sales@cncomon.com)



